



Urine casts: brightfield

Seven cast patterns, identifying features, and common look-alikes

Casts form in renal tubules in a uromodulin protein matrix. Identify the matrix first; loose cell clumps can mimic casts.

Brightfield: transparent casts have little contrast. Reduced illumination helps reveal their outline; phase contrast can clarify cell detail.

Photo	Cast and identifying feature	Association and look-alike
	Hyaline · 400× Pale matrix with parallel sides; faint edges can disappear in strong illumination.	Concentrated urine, exercise, or renal disease; a few may occur without disease. Look-alike: Mucus tapers and varies in width.
	Red cell · 400× Red cells held within a continuous matrix.	Bleeding within the nephron, especially from a glomerular source. Look-alike: Free red cell clumps have no enclosing cast matrix.
	White cell · 400× Granular white cells in a matrix. Nuclear lobes may be difficult to resolve.	Renal inflammation, including pyelonephritis and interstitial nephritis. Look-alike: Tubular cells or white cell clumps can resemble a cast.
	Epithelial / white cell · 400× The source labels this as an epithelial/white-cell cast. Indistinct nuclei limit cell typing.	Tubular injury or renal inflammation; the cast alone cannot establish the cause. Look-alike: White cell and tubular cell casts overlap; compare nuclei, not size alone.
	Granular · 400× Coarse or fine granules fill the matrix as cells and proteins break down.	Numerous muddy brown casts support tubular injury. A few granular casts are less specific. Look-alike: Debris or crystal clumps lack a continuous cast outline.
	Fatty · 400× Refractile lipid droplets within a matrix; some show Maltese crosses in polarized light.	Lipiduria, often with heavy proteinuria. Look-alike: Starch also forms crosses; confirm droplets within a cast.
	Waxy · 400× Highly refractile, sharply outlined matrix; cracks, notches, and blunt ends.	Prolonged tubular stasis, often with advanced renal impairment. Look-alike: Fibers have irregular edges and lack a cast matrix.

Scan low power; identify contents at higher power. Delay and dilute alkaline urine can destroy casts; centrifugation can break or lose them. Broad describes width, not contents. Interpret with other urine findings using the laboratory's reporting method.

Photos: Bert Grijzen, *Microscopic Urine Examination* (2026). CC BY-SA 4.0. Cropped; mixed cast rotated; colors unchanged. Magnifications are source values, not print scale.

Full explanation and sources: latticeMLS.org/guides/urine-casts-brightfield/

Source review: September 5, 2026 · LatticeMLS